

# PRO FREEZER

## LOW TEMPERATURE FREEZER

(-60 °C)



Rack system and boxes available

### FEATURES:

- SINGLE COMPRESSOR/LOW ENERGY
- LESS SPACE REQUIRED
- LOW PRICE
- DIGITAL CONTROLLER
- ENVIRONMENTAL FRIENDLY REFRIGERANT
- REMOTE ALARM (OPTION)
- TEMPERATURE RECORDER (OPTION)



QUALITY SYSTEM AND ENVIRONMENTAL MANAGEMENT SYSTEM  
DS/EN ISO 9001  
DS/EN ISO 14001



Latest update 22.06.2010



**TECHNICAL SPECIFICATIONS:**

		PRO 11	PRO 21	PRO 31	PRO 41	PRO 60
Ambient temperature	°C	+10° to +30°	+10° to +30°	+10° to +30°	+10° to +30°	+10° to +30°
Net. Volume	Litres	130	230	300	360	485
Height with 50 mm castors	mm	920	920	920	920	920
Height with open lid	mm	1620	1620	1620	1620	1620
Width	mm	720	1050	1300	1500	1700
Insulations Thickness	mm	100	100	100	100	100
Depth excl. Handle and hinges	mm	730	730	730	730	730
Power consumption at + 25°C	Watts	330	410	560	610	610
Fuse	A	10	10	10	10	10
Temperature Range	°C	-30° to -60°	-30° to -60°	-30° to -60°	-30° to -60°	-30° to -60°
Noise level	dBa	< 51	< 51	< 51	< 51	< 51
Weight	kg	53	68	78	89	95